

**MATERIAL SAFETY DATA SHEET
FOR
THORN SMITH LABORATORIES
ANALYZED QUANTITATIVE UNKNOWN**

Cement
Catalog Number 80-1150

Manufacturer: Auric Enterprises, Inc.
d/b/a Thorn Smith Laboratories
Address: 7755 Narrow Gauge Road
Beulah, MI 49617
Phone Number: 231-882-4672
MSDS Number: TSL-017
Date Prepared: February 16, 1994
Date Updated: June 23, 2009

SECTION 1 - MATERIAL IDENTIFICATION AND INFORMATION

COMPONENTS - Chemical Name & Common Names
(Hazardous Components 1% or greater; Carcinogens 0.1% or greater)

Portland Cement

Formula: $(4\text{CaO}\cdot\text{Al}_2\text{O}_3\cdot\text{Fe}_2)_3; (3\text{CaOAl}_2\text{O}_3); (3\text{CaO}\cdot\text{SiO}_2); (2\text{CaOSiO}_2)$
CAS No.: N/A
OSHA PEL: N/A
ACGIH TLV: 30 mppcf or 10 mg/m³

SECTION 2 - PHYSICAL/CHEMICAL CHARACTERISTICS

Boiling Point: N/A
Specific Gravity (H₂O=1): 3.0 - 3.2
Vapor Pressure (mm Hg and Temperature): N/A
Melting Point: N/A
Vapor Density (Air-1): N/A
Evaporation Rate (-1): N/A
Solubility in Water: Not Soluble
Appearance and Odor: Gray powder. Odorless.

SECTION 3 - FIRE AND EXPLOSION HAZARD DATA

Flash Point and Method Used: N/A

Auto-Ignition Temperature: N/A

Flammability Limits in Air (% by volume): N/A

LEL: N/A UEL: N/A

Extinguisher Media: Use extinguishing media appropriate for surrounding fire.

Special Fire Fighting Procedures: None identified.

Unusual Fire and Explosion Hazards: None identified.

SECTION 4 - REACTIVITY HAZARD DATA

STABILITY: Stable Unstable

Conditions to Avoid: Stable under ordinary conditions of use and storage.

Incompatibility (Materials to Avoid): N/A

Hazardous Decomposition Products: None identified.

HAZARDOUS POLYMERIZATION: May Occur Will Not Occur

SECTION 5 - HEALTH HAZARD DATA

PRIMARY ROUTES OF ENTRY:

Inhalation Ingestion Skin Contact Eye Contact Not Hazardous

CARCINOGEN LISTED IN:

NTP OSHA IARC Monograph Not Listed

TOXICITY: N/A

HEALTH HAZARDS - Acute:

Skin Contact: May cause irritation to those with sensitive skin.

Eye Contact: May cause irritation.

HEALTH HAZARDS - Chronic: N/A

Signs and Symptoms of Exposure: No information found.

Medical Conditions Generally Aggravated by Exposure: N/A

Other Precautions: Avoid prolonged contact between skin surfaces and wet or moist Portland cement or clothing saturated with wet or moist Portland cement.

EMERGENCY FIRST AID PROCEDURES - Seek medical assistance for further treatment, observation, and support if necessary.

Inhalation: Remove to fresh air.

Skin Contact: Remove any contaminated clothing. Wash promptly and thoroughly with fresh, clean water.

Eye Contact: In case of eye contact, immediately flush with plenty of water for at least 15 minutes, lifting upper and lower eyelids occasionally. Get medical attention immediately.

The information published in this Material Safety Data Sheet has been compiled from our experience and data presented in various technical publications. It is the user's responsibility to determine the suitability of this information for the adoption of necessary safety precautions. We make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and assume no liability resulting from its use. We reserve the right to revise Material Safety Data Sheets periodically as new information becomes available.