

**MATERIAL SAFETY DATA SHEET
FOR
THORN SMITH LABORATORIES
ANALYZED QUANTITATIVE UNKNOWNNS**

Nickel (German) Silver
Catalog Number 80-1385

Manufacturer: Auric Enterprises, Inc.
d/b/a Thorn Smith Laboratories
Address: 7755 Narrow Gauge Road
Beulah, MI 49617
Phone Number: 231-882-4672
MSDS Number: TSL-025
Date Prepared: December 12, 1980
Date Updated: June 23, 2009

SECTION 1 - MATERIAL IDENTIFICATION AND INFORMATION

COMPONENTS - Chemical Name & Common Name
(Hazardous Components 1% or greater; Carcinogens 0.1% or greater)

Copper

Formula: Cu
CAS No.: 7440-50-8
Formula Weight: 63.54
Synonyms: Bronze powder, copper bronze
OSHA PEL: 1 mg/m³ (mist and dust); 0.01 mg/m³ (fume)
ACGIH TLV: 1 mg/m³ (mist and dust); 0.02 mg/m³ (fume)
OTHER LIMITS: N/A

Zinc

Formula: Zn
CAS No.: 7440-66-6
OSHA PEL: N/E
ACGIH TLV: 5 mg/m³
OTHER LIMITS: N/A

Nickel

Formula: Ni
CAS No.: 7440-02-0
OSHA PEL: 1 mg/m³
ACGIH TLV: 1 mg/m³
OTHER LIMITS: N/A

SECTION 2 - PHYSICAL/CHEMICAL CHARACTERISTICS

Boiling Point: N/A

Specific Gravity (H₂O=1): N/A

Vapor Pressure (mm Hg and Temperature): 1 mm @ 1628° (Cu)

Melting Point: 1083° C (Cu)

Vapor Density (Air-1): 8.92 (Cu)

Evaporation Rate (-1): N/A

Solubility in Water: N/A

Water Reactive: No

Appearance and Odor: Light to dark gray color. Fine powder. Distinct reddish metal (Cu).

SECTION 3 - FIRE AND EXPLOSION HAZARD DATA

Flash Point and Method Used: Metal products are not a fire hazard. However, dust generated in grinding operations may present a fire or explosion hazard.

Auto-Ignition Temperature: N/A

Flammability Limits in Air (% by Volume): N/A

LEL: N/A UEL: N/A

Extinguisher Media: Special mixtures of dry chemical suitable for metal fires.

Special Fire Fighting Procedures: Firefighters should wear proper protective equipment and self-contained breathing apparatus with full facepiece operated in the positive pressure mode.

Unusual Fire and Explosion Hazards: N/A

SECTION 4 - REACTIVITY HAZARD DATA

STABILITY: Stable Unstable

Conditions to Avoid: No information found.

Incompatibility (Materials to Avoid): Welding and brazing fumes and gases cannot be classified simply.

Hazardous Decomposition Products: N/A

HAZARDOUS POLYMERIZATION: May Occur Will Not Occur

SECTION 5 - HEALTH HAZARD DATA

PRIMARY ROUTES OF ENTRY:

Inhalation Ingestion Skin Contact Eye Contact Not Hazardous

CARCINOGEN LISTED IN:

NTP OSHA IARC Monograph Not Listed

TOXICITY: Poison to humans by ingestion. (Cu)

HEALTH HAZARDS - Acute:

Inhalation: Excessive inhalation of fumes, gases or dust from many metals can produce an acute reaction known as "metal fume fever." Symptoms consist of chills and fever (very similar to and easily confused with flu symptoms) which come on a few hours after large exposures. May cause headache, nausea and irritation of the upper respiratory tract and lungs. Long-term effects of metal fume fever have not been noted. May cause pulmonary irritation and pneumoconiosis.

Ingestion: Poison to humans by ingestion. (Cu)

Skin Contact: May cause irritation.

Eye Contact: N/A

HEALTH HAZARDS - Chronic:

Inhalation: Industrial exposure to copper fumes, dusts or mists results in metal fume fever with atrophic changes in nasal mucous membranes.

Skin Contact: Certain Nickel compounds have been implicated to be carcinogenic.

Eye Contact: N/A

Ingestion: Metallic taste in mouth, dryness or irritations of the throat followed by coughing, shortness of breath, nausea, fever, body ache, and chills.

Signs and Symptoms of Exposure: Short term symptoms may include a metallic taste in the mouth, dryness or irritations of the throat followed by coughing, shortness of breath, nausea, fever, body ache, and chills. Inhalation of fumes may cause headache, nausea, and irritation of the nose, throat, and lungs. Individuals with impaired pulmonary functions or illness may have symptoms exacerbated by fume irritants.

Medical Conditions Generally Aggravated by Exposure: See Above.

EMERGENCY FIRST AID PROCEDURES - Seek medical assistance for further treatment, observation, and support if necessary.

Inhalation: If acute overexposure to fumes occurs, move the exposed person from the adverse environment and seek medical attention.

Ingestion: If particles are ingested, give 1 to 2 glasses of water or milk. Induce vomiting only if victim is fully conscious and has not convulsed. All ingestion cases should have immediate medical aid.

Skin Contact: Remove any contaminated clothing. Wipe off excess from skin. Immediately wash skin with soap and water for at least 15 minutes. Get medical attention if irritation develops or persists.

Eye Contact: Immediately flush with plenty of water for at least 15 minutes, lifting upper and lower eyelids occasionally.

SARA/TITLE III HAZARD CATEGORIES AND LISTS

Acute: Yes Chronic: Yes Flammability: Yes Pressure: No Reactivity: No

Extremely Hazardous Substance: No

CERCLA Hazardous Substance: Yes

SARA 313 Toxic Chemicals: Yes

TSCA Inventory: No (Cu)

SECTION 6 - CONTROL AND PROTECTIVE MEASURES

Respiratory Protection (Specify Type): Use appropriate NIOSH approved respirators where local exhaust or ventilation does not keep exposure below TLV.

Protective Gloves: As needed to protect against physical hazards.

Eye Protection: Wear chemical safety goggles.

VENTILATION TO BE USED: General or local exhaust ventilation to keep fume or dust levels as low as possible.

Local Exhaust Mechanical (General) Special

Other (Specify)

Other Protective Clothing and Equipment: Wear clean body-covering clothing.

Hygienic Work Practices: Avoid contact with eyes, skin, and clothing. Avoid breathing dust. Keep container closed when not in use. Use with adequate ventilation. Wash thoroughly after handling. Keep away from incompatibles.

SECTION 7 - PRECAUTIONS FOR SAFE HANDLING & USE/LEAK PROCEDURES

Steps to be taken if material is spilled or released: N/A

Waste Disposal Methods: Dispose in accordance with all applicable local, state, and federal environmental regulations.

Precautions to be taken in handling and storage: Keep in a tightly closed container. Suitable for any general chemical storage area. Store in accordance with all applicable local, state, and federal environmental regulations.

Other precautions and/or special hazards: No information available.

NFPA Rating: Mixture not rated.

HMIS Rating: No information available.

SECTION 8 - TRANSPORTATION DATA AND ADDITIONAL INFORMATION

Domestic (D.O.T.)

Proper Shipping Name: Chemicals, n.o.s. (Non-regulated)

International (T.M.O.)

Proper Shipping Name: Chemicals, n.o.s. (Non-regulated)

AIR (I.C.A.O.)

Proper Shipping Name: Chemicals, n.o.s. (Non-regulated)

The information published in this Material Safety Data Sheet has been compiled from our experience and data presented in various technical publications. It is the user's responsibility to determine the suitability of this information for the adoption of necessary safety precautions. We make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and assume no liability resulting from its use. We reserve the right to revise Material Safety Data Sheets periodically as more information becomes available.